



**HEALTH GORILLA™**  
Secure Clinical Network

# Health Gorilla Clinical Network™ API

## HEALTH GORILLA API LICENSE AGREEMENT

In order to use the Health Gorilla Clinical Network API a legal agreement must be accepted. The legal agreement (the “Agreement”) between the you and, if applicable, your company, organization or other legal entity for whom you have authority to enter into this Agreement (“Developer”) and Health Gorilla, Inc. (“Health Gorilla”), for use of certain of Health Gorilla's Application Programming Interfaces (“APIs”) that enable Developer to provide additional functionality for users of Health Gorilla products or to integrate access to certain Health Gorilla products into Developer's product offerings. By clicking “I agree to the terms of service,” or otherwise accessing the APIs, Developer indicates its acceptance of this Agreement.

Last revised September 21 2018

# Table of Contents

<b>Health Gorilla Clinical Network™ API</b>	<b>1</b>
<b>Table of Contents</b>	<b>2</b>
<b>FHIR RESTful Health Gorilla API</b>	<b>4</b>
<b>Currently supported Endpoints</b>	<b>4</b>
<b>Organization</b>	<b>5</b>
Endpoint	5
FHIR Operations	5
Examples	5
<b>Practitioner</b>	<b>6</b>
Endpoint	6
FHIR Operations	6
Examples	6
<b>Patient</b>	<b>7</b>
Endpoint	7
FHIR Operations	7
Examples	8
<b>AllergyIntolerance</b>	<b>9</b>
Endpoint	9
FHIR Operations	9
Examples	9
<b>CarePlan</b>	<b>10</b>
Endpoint	10
FHIR Operations	10
Examples	10
<b>Condition</b>	<b>11</b>
Endpoint	11
FHIR Operations	11
Examples	12
<b>DeviceUseStatement</b>	<b>13</b>
Endpoint	13
FHIR Operations	13
Examples	13
<b>DiagnosticReport</b>	<b>14</b>
Endpoint	14
FHIR Operations	14
Examples	14
<b>DocumentReference</b>	<b>16</b>
Endpoint	16
FHIR Operations	16
Examples	16
<b>Encounter</b>	<b>17</b>
Endpoint	17
FHIR Operations	17
Examples	17

<b>FamilyMemberHistory</b>	<b>18</b>
Endpoint	18
FHIR Operations	18
Examples	18
<b>Immunization</b>	<b>19</b>
Endpoint	19
FHIR Operations	19
Examples	19
<b>MedicationStatement</b>	<b>20</b>
Endpoint	20
FHIR Operations	20
Examples	20
<b>Observation</b>	<b>21</b>
Endpoint	21
FHIR Operations	21
Examples	21
<b>Procedure</b>	<b>22</b>
Endpoint	22
FHIR Operations	22
Examples	22
<b>ProcedureRequest</b>	<b>23</b>
Endpoint	23
FHIR Operations	23
Examples	23
<b>RequestGroup</b>	<b>24</b>
Endpoint	24
FHIR Operations	24
Examples	24
<b>Subscriptions</b>	<b>25</b>
Endpoint	25
FHIR Operations	25
Examples	26
Subscription channel types	26

# FHIR RESTful Health Gorilla API

FHIR®, or Fast Healthcare Interoperability Resources, provides a lightweight REST-based access layer for standard HL7-defined data models.

The most actual list of endpoints can be found at the following URL:

<https://api.healthgorilla.com/fhir/metadata>

## Currently supported Endpoints

- [Organization](#)
- [Practitioner](#)
- [Patient](#)
- [AllergyIntolerance](#)
- [CarePlan](#)
- [Condition](#)
- [DeviceUseStatement](#)
- [DiagnosticReport](#)
- [DocumentReference](#)
- [Encounter](#)
- [FamilyMemberHistory](#)
- [Immunization](#)
- [MedicationStatement](#)
- [Observation](#)
- [Procedure](#)
- [ProcedureRequest](#)
- [RequestGroup](#)
- [Subscriptions](#)

# Organization

Using Organization resource, you can obtain information about Practices, Laboratories, Hospitals and other facilities managed by HealthGorilla.

FHIR Resource: Organization <https://www.hl7.org/fhir/organization.html>

The following attributes are supported:

Field Name	Comments
<code>active</code>	
<code>type</code>	There are several value sets supported: <ol style="list-style-type: none"><li>1. <a href="https://www.hl7.org/fhir/valueset-organization-type.html">https://www.hl7.org/fhir/valueset-organization-type.html</a></li><li>2. <a href="https://www.healthgorilla.com/fhir/ValueSet/facility-type">https://www.healthgorilla.com/fhir/ValueSet/facility-type</a></li><li>3. <a href="https://www.healthgorilla.com/fhir/ValueSet/organization-type">https://www.healthgorilla.com/fhir/ValueSet/organization-type</a></li></ol>
<code>name</code>	
<code>telecom</code>	
<code>address</code>	
<code>partOf</code>	Link to the department or to the parent Organization.

## Endpoint

<https://api.healthgorilla.com/fhir/Organization>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/Organization/{ID}</code>	GET	Retrieve organization by ID

### 2. Search

HTTP Request	Method	Action
<code>/Organization?[parameter={value}]</code>	GET	Retrieve Organizations by the specified search criteria

### Search Parameters:

Name	Type	Description
type	Token	A code for the type of organization.
name	string	A portion of the organization's name or alias.
address-state	string	A state specified in an address.
address-city	string	A city specified in an address.
address-postalcode	string	A postal code specified in an address.
ordering-enabled	token	A ordering-enabled property. Use to get only organizations that support ordering.
_count	number	A positive integer. If passed, then server will not return more than specified.

### Examples:

1. Search laboratories with name LabCorp that support ordering:

```
https://api.healthgorilla.com/fhir/Organization?type=http://hl7.org/fhir/organization-type|prov&name=LabCorp&ordering-enabled=true
```

2. Search practices in Sunnyvale:

```
https://api.healthgorilla.com/fhir/Organization?type=https://www.healthgorilla.com/fhir/organization-type|PR&address-city=Sunnyvale&name=HealthCenter
```

3. Search AETNA insurance companies in California

```
https://api.healthgorilla.com/fhir/Organization?type=http://hl7.org/fhir/organization-type|ins&name=AETNA&address-state=CA&_count=15
```

## Practitioner

Using Practitioner resource, you can obtain information about Healthcare providers managed by HealthGorilla.

FHIR Resource: Practitioner <https://www.hl7.org/fhir/practitioner.html>

The following attributes are supported:

Field Name	Comments
------------	----------

identifier	- NPI - UID
active	
birthDate	
name	
telecom	
address	
qualification	
<a href="http://hl7.org/fhir/extension-practitioner-classification.html">http://hl7.org/fhir/extension-practitioner-classification.html</a>	Administrative Role
<a href="http://hl7.org/fhir/extension-event-onbehalfof.html">http://hl7.org/fhir/extension-event-onbehalfof.html</a>	Parent Organization

## Endpoint

<https://api.healthgorilla.com/fhir/Practitioner>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/Practitioner/{ID}</code>	GET	Retrieve provider by ID

### 2. Search

HTTP Request	Method	Action
<code>/Practitioner?[parameter={value}]</code>	GET	Retrieve providers by the specified search criteria

Search Parameters:

Name	Type	Description
Identifier	Token	A practitioner's Identifier.

Examples:

1. Search practitioner by NPI number:

```
https://api.healthgorilla.com/fhir/Practitioner?identifier=http://hl7.org/fhir/sid/us-npi|1234567893
```

## 2. Search practitioner by HealthGorilla assigned identifier:

```
https://api.healthgorilla.com/fhir/Practitioner?identifier=https://www.healthgorilla.com|f3c88f5b53910aba83b224f4
```

## Patient

The FHIR Patient resource defines demographics, care providers, and other administrative information about a person receiving care at a health organization. The Patient resource can be thought of as the starting point for many different client application workflows; often, a workflow will begin with a patient query and selection, and the server ID will be used to drive additional queries for diagnostic and care-related information.

FHIR Resource: Patient <https://www.hl7.org/fhir/patient.html>

## Endpoint

<https://api.healthgorilla.com/fhir/Patient>

The following attributes are supported:

Field Name	Comments
active	
identifier	- MRN - SSN
photo	
name	
gender	
birthDate	
telecom	
address	
maritalStatus	
communication	
managingOrganization	



link	
<a href="http://hl7.org/fhir/extension-patient-interpreterrequired.html">http://hl7.org/fhir/extension-patient-interpreterrequired.html</a>	Interpreter Required
<a href="http://hl7.org/fhir/us/core/StructureDefinition-us-core-race.html">http://hl7.org/fhir/us/core/StructureDefinition-us-core-race.html</a>	Race
<a href="http://hl7.org/fhir/us/core/StructureDefinition-us-core-ethnicity">http://hl7.org/fhir/us/core/StructureDefinition-us-core-ethnicity</a>	Ethnicity

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/Patient/{ID}</code>	GET	Retrieve patient by ID

### 2. Search

HTTP Request	Method	Action
<code>/Patient?[parameter={value}]</code>	GET	Retrieve Patient resources meeting the specified search criteria

Search Parameters:

Name	Type	Description
identifier	Token	Patient Identifier (MRN)
family	String	Patient Family Name
given	String	Patient Given Name
birthdate	String	Patient DOB (in YYYYMMDD format)
gender	Token	Patient Gender (F/M)
address-postalcode	String	Patient Postal (zip) code

Examples:

#### 1. Search patient by MRN:

```
https://api.healthgorilla.com/fhir/Patient?identifier=https://www.healthgorilla.com|8063965523692105bfc35def
```

#### 2. Search patient by SSN:

```
https://api.healthgorilla.com/fhir/Patient?identifier=http://hl7.org/fhir/sid/us-ssn|000-00-0001
```

### 3. Search patient by name:

```
https://api.healthgorilla.com/fhir/Patient?given=Homer&family=Simpson
```

### 4. Search males were born after October 1998 and whose firstname contains 'Alex':

```
https://api.healthgorilla.com/fhir/Patient?given:contains=Alex&birthdate=gt1998-10&gender=male
```

### 3. Create

HTTP Request	Method	Action
/Patient	POST	Creates patient resource

The following attributes are required:

Name	Details
name	
gender	
birthDate	

#### 3.1. Conditional create

This method also supports [conditional create](#). To use it, you need to set additional header:

If-None-Exist: [search parameters]

Supported parameters:

- given
- family
- gender
- birthdate

Example:

If-None-Exist:

```
given=Rajesh&family=Koothrappali&gender=male&birthdate=1981-01-30
```

### 4. Update

HTTP Request	Method	Action
/Patient/{ID}	PUT	Updates patient resource

# AllergyIntolerance

The FHIR AllergyIntolerance resource defines clinical information about a patient's allergic response to a substance. The AllergyIntolerance resource defines the substance that elicited the response, as well as when the reaction occurred, the severity, and the type of reaction noted.

FHIR Resource: AllergyIntolerance <https://www.hl7.org/fhir/allergyintolerance.html>

The following attributes are supported:

Name	Comment
<a href="#">verificationStatus</a>	
<a href="#">criticality</a>	
<a href="#">type</a>	
<a href="#">category</a>	
<a href="#">clinicalStatus</a>	
<a href="#">code</a>	
<a href="#">patient</a>	
<a href="#">recorder</a>	
<a href="#">onset</a>	
<a href="#">assertedDate</a>	Documentation Date
<a href="#">lastOccurrence</a>	Completion Date
<a href="#">reaction</a>	
<a href="#">reaction.manifestation</a>	
<a href="#">reaction.severity</a>	
<a href="#">reaction.substance</a>	
<a href="#">reaction.description</a>	

## Endpoint

<https://api.healthgorilla.com/fhir/AllergyIntolerance>

## FHIR Operations

The following operations are currently supported:

## 1. Read

HTTP Request	Method	Action
<code>/AllergyIntolerance/{ID}</code>	GET	Retrieve allergy resource by ID

## 2. Search

HTTP Request	Method	Action
<code>/AllergyIntolerance?[parameter={value}]</code>	GET	Retrieve allergies by the specified search criteria

Search Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

1. Retrieve allergies by patient:

```
https://api.healthgorilla.com/fhir/AllergyIntolerance?patient=8063965523692105bfc35def
```

## 3. Create

HTTP Request	Method	Action
<code>/AllergyIntolerance</code>	POST	Creates allergy resource

The following attributes are required:

Name	Details
<code>type</code>	
<code>category</code>	
<code>code</code>	
<code>patient</code>	

## 4. Update

HTTP Request	Method	Action
<code>/AllergyIntolerance/{ID}</code>	PUT	Updates allergy resource or creates new one

The following attributes are immutable:

Name	Details
------	---------

type	
category	
patient	
code	
recorder	

## 5. Delete

HTTP Request	Method	Action
<code>/AllergyIntolerance/{ID}</code>	DELETE	Deletes allergy resource

# CarePlan

The CarePlan resource describes the intention of how one or more practitioners intend to deliver care for a patient, group or community for a period of time, possibly limited to care for a specific condition or set of conditions.

FHIR Resource: CarePlan <https://www.hl7.org/fhir/careplan.html>

The following attributes are supported:

Name	Comment
subject	
status	
description	
participant	
careTeam	
period	
intent	
goal	
activity	
activity.outcomeReference	MedicationStatement, Encounter, Procedure, DeviceUseStatement
activity.detail	
activity.detail.category	
activity.detail.code	
activity.detail.description	
activity.detail.scheduled	
activity.detail.quantity	
activity.detail.product	
activity.detail.location	
activity.detail.performer	

activity.detail.reasonCode	
activity.detail.reasonReference	
note	
<a href="https://www.healthgorilla.com//fhir/StructureDefinition/record-recorder">https://www.healthgorilla.com//fhir/StructureDefinition/record-recorder</a>	Recorder

## Endpoint

<https://api.healthgorilla.com/fhir/CarePlan>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<a href="#">/CarePlan/{ID}</a>	GET	Retrieve CarePlan resource by ID

### 2. Search

HTTP Request	Method	Action
<a href="#">/CarePlan?[parameter={value}]</a>	GET	Retrieve CarePlan resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

#### 1. Retrieve CarePlans by patient:

```
https://api.healthgorilla.com/fhir/CarePlan?patient=8063965523692105bfc35def
```

### 3. Create

HTTP Request	Method	Action
<a href="#">/CarePlan</a>	POST	Creates CarePlan resource

Activity may be passed using [activity.outcomeReference](#) or inside [activity.detail](#).

The following attributes are required:

Name	Details
subject	

#### 4. Update

HTTP Request	Method	Action
/CarePlan/{ID}	PUT	Updates CarePlan resource or creates new one

The following attributes are immutable:

Name	Details
subject	
recorder	

#### 5. Delete

HTTP Request	Method	Action
/CarePlan/{ID}	DELETE	Deletes CarePlan resource



## Condition

The Condition resource describes a clinical condition, problem, diagnosis, or other event, situation, issue, or clinical concept that has risen to a level of concern.

FHIR Resource: Condition <https://www.hl7.org/fhir/condition.html>

The following attributes are supported:

Name	Comment
<a href="#">clinicalStatus</a>	
<a href="#">onset</a>	
<a href="#">subject</a>	
<a href="#">code</a>	ICD10, Snomed
<a href="#">asserter</a>	
<a href="#">assertedDate</a>	

## Endpoint

<https://api.healthgorilla.com/fhir/Condition>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/Condition/{ID}</code>	GET	Retrieve Condition resource by ID

### 2. Search

HTTP Request	Method	Action
<code>/Condition?[parameter={value}]</code>	GET	Retrieve Condition resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

1. Retrieve Conditions by patient:

```
https://api.healthgorilla.com/fhir/Condition?patient=8063965523692105bfc35def
```

### 3. Create

HTTP Request	Method	Action
<code>/Condition</code>	POST	Creates Condition resource

The following attributes are required:

Name	Details
<code>subject</code>	
<code>code</code>	

### 4. Update

HTTP Request	Method	Action
<code>/Condition/{ID}</code>	PUT	Updates Condition resource

The following attributes are immutable:

Name	Details
<code>subject</code>	
<code>code</code>	

### 5. Delete

HTTP Request	Method	Action
<code>/Condition/{ID}</code>	DELETE	Deletes Condition resource

## Coverage

The Coverage resource describes a patient's insurance.

FHIR Resource: Condition <https://www.hl7.org/fhir/coverage.html>

The following attributes are supported:

Name	Comment
identifier	Policy ID
type	
beneficiary	Related patient
relationship	Self, Spouse, Parent, Other
subscriber	
payor	Reference to insurance provider
grouping	
grouping.group	Group Number
grouping.planDisplay	Plan Name
order	1 - primary insurance 2 - secondary insurance 3 - tertiary insurance

## Endpoint

<https://api.healthgorilla.com/fhir/Coverage>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<a href="#">/Coverage/{ID}</a>	GET	Retrieve Coverage resource by ID

### 2. Search

HTTP Request	Method	Action
--------------	--------	--------

<code>/Coverage?[parameter={value}]</code>	GET	Retrieve Condition resources by the specified search criteria
--	-----	---

Parameters:

Key	Type	Description
beneficiary	Reference	Patient Identifier

Examples:

1. Retrieve insurances by patient:

```
https://api.healthgorilla.com/fhir/Coverage?beneficiary=8063965523692105bfc35def
```

### 3. Create

HTTP Request	Method	Action
<code>/Coverage</code>	POST	Adds insurance to patient

The following attributes are required:

Name	Details
identifier	
beneficiary	
relationship	
payor	
grouping	

### 4. Delete

HTTP Request	Method	Action
<code>/Coverage/{ID}</code>	DELETE	Deletes insurance from patient

## DeviceUseStatement

The DeviceUseStatement resource is a record of a device being used by a patient.

FHIR Resource: DeviceUseStatement <https://www.hl7.org/fhir/deviceusestatement.html>

The following attributes are supported:

Name	Comment
subject	Patient
whenUsed	
source	Practitioner
indication	
note	
device	

Device Resource support the following attributes:

Name	Comment
identifier	Serial Number (SNO)
type	
manufacturer	
manufactreDate	
model	
expirationDate	
lotNumber	
note	

## Endpoint

<https://api.healthgorilla.com/fhir/DeviceUseStatement>

## FHIR Operations

The following operations are currently supported:

## 1. Read

HTTP Request	Method	Action
<code>/DeviceUseStatement/{ID}</code>	GET	Retrieve DeviceUseStatement resource by ID

## 2. Search

HTTP Request	Method	Action
<code>/DeviceUseStatement?[parameter={value}]</code>	GET	Retrieve DeviceUseStatement resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

### 1. Retrieve Conditions by patient:

```
https://api.healthgorilla.com/fhir/DeviceUseStatement?patient=8063965523692105bfc35def
```

## 3. Create

HTTP Request	Method	Action
<code>/DeviceUseStatement</code>	POST	Creates DeviceUseStatement resource

The following attributes are required:

Name	Details
<code>subject</code>	
<code>device</code>	
<code>device.type</code>	

## 4. Update

HTTP Request	Method	Action
<code>/DeviceUseStatement/{ID}</code>	PUT	Updates DeviceUseStatement resource

The following attributes are immutable:

Name	Details
------	---------

subject	
device	

## 5. Delete

HTTP Request	Method	Action
<code>/DeviceUseStatement/{ID}</code>	DELETE	Deletes DeviceUseStatement resource

## DiagnosticReport

The DiagnosticReport resource contains the findings and interpretation of diagnostic tests performed on patients, groups of patients, devices, and locations, and/or specimens derived from these. The report includes clinical context such as requesting and provider information, and some mix of atomic results, images, textual and coded interpretations, and formatted representation of diagnostic reports.

FHIR Resource: DiagnosticReport <https://www.hl7.org/fhir/diagnosticreport.html>

The following attributes are supported:

Name	Comment
identifier	HG ID, ACSN
status	
code	
category	
subject	Patient
effective	
issued	
performer	Organization
result	Observations
conclusion	
codedDiagnosis	
presentedForm	

basedOn	ProcedureRequests
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/patient-location">https://www.healthgorilla.com/fhir/StructureDefinition/patient-location</a>	Patient Location (Room, Bad, Hospital)

## Endpoint

<https://api.healthgorilla.com/fhir/DiagnosticReport>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<a href="#">/DiagnosticReport/{ID}</a>	GET	Retrieve DiagnosticReport resource by ID

### 2. Search

HTTP Request	Method	Action
<a href="#">/DiagnosticReport?[parameter={value}]</a>	GET	Retrieve DiagnosticReport resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier
category	Token	Which diagnostic discipline/department created the report
subject	Reference	The subject of the report
based-on	Reference	Reference to the procedure request.
_lastUpdated	Date	

Examples:

#### 1. Retrieve diagnostic reports by patient:

```
https://api.healthgorilla.com/fhir/DiagnosticReport?patient=8063965523692105bfc35def
```

#### 2. Get all laboratory reports by updated since 2017:

```
https://api.healthgorilla.com/fhir/DiagnosticReport?category=http://hl7.org/fhir/v2/0074|LAB&_lastUpdated=gt2017
```



## DocumentReference

A DocumentReference resource is used to describe a document that is made available to a healthcare system. See: <https://www.hl7.org/fhir/documentreference.html>

FHIR Resource: DocumentReference <https://www.hl7.org/fhir/documentreference.html>

The following attributes are supported:

Name	Comment
<a href="#">identifier</a>	
<a href="#">docStatus</a>	
<a href="#">type</a>	<ul style="list-style-type: none"><li>- Loinc</li><li>- <a href="https://www.healthgorilla.com/document-type">https://www.healthgorilla.com/document-type</a></li></ul>
<a href="#">class</a>	
<a href="#">subject</a>	Patient
<a href="#">created</a>	
<a href="#">description</a>	
<a href="#">indexed</a>	
<a href="#">author</a>	Organization
<a href="#">authenticator</a>	Practitioner
<a href="#">custodian</a>	
<a href="#">securityLabel</a>	
<a href="#">content</a>	
<a href="#">content.attachment</a>	
<a href="#">related</a>	

## Endpoint

<https://api.healthgorilla.com/fhir/DocumentReference>

## FHIR Operations

The following operations are currently supported:

## 1. Read

HTTP Request	Method	Action
<code>/DocumentReference/{ID}</code>	GET	Retrieve document by ID

## 2. Search

HTTP Request	Method	Action
<code>/DocumentReference?[parameter={value}]</code>	GET	Retrieve documents by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier
type	Token	Kind of document (LOINC if possible)
class	Token	Categorization of document
_lastUpdated	Date	

Examples:

1. Retrieve documents by patient:

```
https://api.healthgorilla.com/fhir/DocumentReference?patient=8063965523692105bfc35def
```

2. Get all Discharge Summary documents updated in August 2018:

```
https://api.healthgorilla.com/fhir/DocumentReference?type=18842-5&_lastUpdated=2018-08
```

## 3. Create

HTTP Request	Method	Action
<code>/DocumentReference</code>	POST	Creates new Document to the specified patient chart

The following attributes are required:

Name	Details
<code>subject</code>	
<code>type</code>	
<code>content.attachment</code>	

## Example

```
{
  "resourceType": "DocumentReference",
  "status": "current",
  "type": {
    "coding": [
      {
        "system": "http://loinc.org",
        "code": "56447-6",
        "display": "Plan of care note"
      },
      {
        "system":
"https://www.healthgorilla.com/document-type",
        "code": "PlanOfCare",
        "display": "Plan of Care"
      }
    ],
    "text": "Plan of Care"
  },
  "class": {
    "coding": [
      {
        "system": "http://loinc.org",
        "code": "56447-6",
        "display": "Plan of care note"
      }
    ]
  },
  "subject": {
    "reference": "Patient/03db435262ce2aba13b8d8c6",
    "display": "Lisa Simpson"
  },
  "created": "2018-10-02T04:30:00+03:00",
  "author": [
    {
      "display": "El Camino Hospital"
    }
  ],
  "authenticator": {
    "reference": "Practitioner/fc076a524383b188ffc10051",
    "display": "Robert Kelso, M.D."
  },
  "content": [
    {
      "attachment": {
```

```

        "contentType": "text/plain",
        "data": "d2Vla2x5IGRpZ2VzdCBj...",
        "size": 186,
        "title": "file.txt"
    }
}
]
}

```

#### 4. Delete

HTTP Request	Method	Action
/DocumentReference/{ID}	DELETE	Deletes Document

### Encounter

The Encounter resource describes an interaction between a patient and healthcare provider(s) for the purpose of providing healthcare service(s) or assessing the health status of a patient.

FHIR Resource: Encounter <https://www.hl7.org/fhir/encounter.html>

The following attributes are supported:

Name	Comment
type	
subject	Patient
period	
reason	
diagnosis	
participant	
location	
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder">https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder</a>	Recorder (Participant)

## Endpoint

<https://api.healthgorilla.com/fhir/Encounter>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/Encounter/{ID}</code>	GET	Retrieve Encounter resource by ID

### 2. Search

HTTP Request	Method	Action
<code>/Encounter?[parameter={value}]</code>	GET	Retrieve Encounter resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

#### 1. Retrieve encounters by patient:

```
https://api.healthgorilla.com/fhir/Encounter?patient=8063965523692105bfc35def
```

### 3. Create

HTTP Request	Method	Action
<code>/Encounter</code>	POST	Creates Encounter resource

The following attributes are required:

Name	Details
<code>subject</code>	
<code>type</code>	

### 4. Update

HTTP Request	Method	Action
<code>/Encounter/{ID}</code>	PUT	Updates Encounter resource

The following attributes are immutable:

Name	Details
subject	
type	
recorder	

## 5. Delete

HTTP Request	Method	Action
<code>/Encounter/{ID}</code>	DELETE	Deletes Encounter resource

# FamilyMemberHistory

The FamilyMemberHistory resource describes significant health events and conditions for a person related to the patient relevant in the context of care for the patient. This resource returns the current status, age, and medical and surgical history of relatives related to a specific patient.

FHIR Resource: FamilyMemberHistory <https://www.hl7.org/fhir/familymemberhistory.html>

The following attributes are supported:

Name	Comment
<a href="#">patient</a>	
<a href="#">name</a>	
<a href="#">relationship</a>	
<a href="#">gender</a>	
<a href="#">born</a>	
<a href="#">deceased</a>	
<a href="#">condition</a>	
<a href="#">condition.code</a>	
<a href="#">condition.outcome</a>	
<a href="#">condition.onset</a>	
<a href="#">condition.note</a>	
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder">https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder</a>	Recorder (Participant)

## Endpoint

<https://api.healthgorilla.com/fhir/FamilyMemberHistory>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<a href="#">/FamilyMemberHistory/{ID}</a>	GET	Retrieve FamilyMemberHistory resource by ID

## 2. Search

HTTP Request	Method	Action
<code>/FamilyMemberHistory?[parameter={value}]</code>	GET	Retrieve FamilyMemberHistory resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

1. Retrieve FamilyMemberHistory records by patient:

```
https://api.healthgorilla.com/fhir/FamilyMemberHistory?patient=8063965523692105bfc35def
```

## 3. Create

HTTP Request	Method	Action
<code>/FamilyMemberHistory</code>	POST	Creates FamilyMemberHistory resource

The following attributes are required:

Name	Details
<code>patient</code>	
<code>relationship</code>	

## 4. Update

HTTP Request	Method	Action
<code>/FamilyMemberHistory/{ID}</code>	PUT	Updates FamilyMemberHistory resource

The following attributes are immutable:

Name	Details
<code>patient</code>	
<code>recorder</code>	

## 5. Delete



HTTP Request	Method	Action
<code>/FamilyMemberHistory/{ID}</code>	DELETE	Deletes FamilyMemberHistory resource

## Immunization

The Immunization resource describes the event of a patient being administered a vaccination or a record of a vaccination as reported by a patient, a clinician or another party and may include vaccine reaction information and what vaccination protocol was followed.

FHIR Resource: Immunization <https://www.hl7.org/fhir/immunization.html>

The following attributes are supported:

Name	Comment
<code>status</code>	
<code>notGiven</code>	
<code>vaccine</code>	
<code>patient</code>	
<code>date</code>	
<code>deceased</code>	
<code>expirationDate</code>	
<code>doseQuantity</code>	
<code>route</code>	
<code>practitioner</code>	
<code>location</code>	
<code>note</code>	

The following Vaccine attributes are supported:

Name	Comment
<code>vaccineCode</code>	
<code>manufacturer</code>	
<code>lotNumber</code>	

## Endpoint

<https://api.healthgorilla.com/fhir/Immunization>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/Immunization/{ID}</code>	GET	Retrieve Immunization by ID

### 2. Search

HTTP Request	Method	Action
<code>/Immunization?[parameter={value}]</code>	GET	Retrieve immunizations by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

#### 1. Retrieve immunizations records by patient:

```
https://api.healthgorilla.com/fhir/Immunization?patient=8063965523692105bfc35def
```

### 3. Create

HTTP Request	Method	Action
<code>/Immunization</code>	POST	Creates immunization resource

The following attributes are required:

Name	Details
patient	
vaccine	

### 4. Update

HTTP Request	Method	Action
--------------	--------	--------

<code>/Immunization/{ID}</code>	PUT	Updates immunization resource
---------------------------------	-----	-------------------------------

The following attributes are immutable:

Name	Details
<code>patient</code>	
<code>practitioner</code>	
<code>vaccine</code>	

## 5. Delete

HTTP Request	Method	Action
<code>/Immunization/{ID}</code>	DELETE	Deletes immunization resource

## MedicationStatement

The MedicationStatement resource describes a record of a medication that is being consumed by a patient. A MedicationStatement may indicate that the patient may be taking the medication now or has taken the medication in the past or will be taking the medication in the future.

FHIR Resource: MedicationStatement <https://www.hl7.org/fhir/medicationstatement.html>

The following attributes are supported:

Name	Comment
<code>status</code>	
<code>effective</code>	
<code>medication</code>	
<code>subject</code>	Patient
<code>dosage</code>	
<code>reasonCode</code>	
<code>reasonReference</code>	
<code>partOf</code>	MedicationDispense
<code>taken</code>	
<code>note</code>	

<a href="https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder">https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder</a>	Recorder (Participant)
---	------------------------

The following Medication attributes are supported:

Name	Comment
code	
manufacturer	
ingredient	
form	

The following Dosage attributes are supported:

Name	Comment
dose	
timing	
route	
additionalInstruction	

The following MedicationDispense attributes are supported:

Name	Comment
medication	
subject	Patient
quantity	
performer	
note	
dosageInstruction	
<a href="http://hl7.org/fhir/StructureDefinition/medicationdispense-validityPeriod">http://hl7.org/fhir/StructureDefinition/medicationdispense-validityPeriod</a>	

## Endpoint

<https://api.healthgorilla.com/fhir/MedicationStatement>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/MedicationStatement/{ID}</code>	GET	Retrieve MedicationStatement resource by ID

### 2. Search

HTTP Request	Method	Action
<code>/MedicationStatement?[parameter={value}]</code>	GET	Retrieve MedicationStatement resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

1. Retrieve MedicationStatement records by patient:

```
https://api.healthgorilla.com/fhir/MedicationStatement?patient=8063965523692105bfc35def
```

### 3. Create

HTTP Request	Method	Action
<code>/MedicationStatement</code>	POST	Creates MedicationStatement resource

The following attributes are required:

Name	Details
<code>subject</code>	
<code>medication</code>	
<code>taken</code>	

### 4. Update

HTTP Request	Method	Action
<code>/MedicationStatement/{ID}</code>	PUT	Updates MedicationStatement resource

The following attributes are immutable:

Name	Details
patient	
recorder	
medication	

## 5. Delete

HTTP Request	Method	Action
<code>/MedicationStatement/{ID}</code>	DELETE	Deletes MedicationStatement resource

## Observation

The Observation resource describes measurements and simple assertions made about a patient, device or other subject.

FHIR Resource: Observation <https://www.hl7.org/fhir/observation.html>

HealthGorilla supports the following categories of Observation:

System	Code	Display
<a href="http://hl7.org/fhir/observation-category">http://hl7.org/fhir/observation-category</a>	laboratory	Laboratory
<a href="http://hl7.org/fhir/observation-category">http://hl7.org/fhir/observation-category</a>	vital-signs	Vital Sign
<a href="http://hl7.org/fhir/observation-category">http://hl7.org/fhir/observation-category</a>	social-history	Social History
<a href="https://www.healthgorilla.com/fhir/observation-category">https://www.healthgorilla.com/fhir/observation-category</a>	functional-status	Functional Status

The following attributes are supported:

Name	Comment
status	
category	
code	
subject	Patient
effective	
issued	
performer	
value	
interpretation	
comment	
referenceRange	

<a href="#">related</a>	Observation, DeviceUseStatement
<a href="http://hl7.org/fhir/StructureDefinition/event-partOf">http://hl7.org/fhir/StructureDefinition/event-partOf</a>	DiagnosticReport, Observation

## Endpoint

<https://api.healthgorilla.com/fhir/Observation>

## FHIR Operations

The following operations are currently supported:

List of available operations per category:

Laboratory	Read, Search
Vital Sign	Read, Search, Create, Update, Delete
Social History	Read, Search, Create, Update, Delete
Functional Status	Read, Search, Create, Update, Delete

### 1. Read

HTTP Request	Method	Action
<a href="#">/Observation/{ID}</a>	GET	Retrieve Observation resource by ID

### 2. Search

HTTP Request	Method	Action
<a href="#">/Observation?[parameter={value}]</a>	GET	Retrieve Observation resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier
category	Token	The classification of the type of observation (see: <a href="http://hl7.org/fhir/codesystem-observation-category.html">http://hl7.org/fhir/codesystem-observation-category.html</a> )

Examples:

1. Retrieve Laboratory Observations by patient:

```
https://api.healthgorilla.com/fhir/Observation?category=http://hl7.org/fhir/observation-category|laboratory
```



```
&patient=8063965523692105bfc35def
```

## 2. Retrieve Vital Signs by patient:

```
https://api.healthgorilla.com/fhir/Observation?  
category=http://hl7.org/fhir/observation-category|vital-signs  
&patient=8063965523692105bfc35def
```

## 3. Retrieve Social History by patient:

```
https://api.healthgorilla.com/fhir/Observation?  
category=http://hl7.org/fhir/observation-category|social-history  
&patient=8063965523692105bfc35def
```

## 4. Retrieve Functional Status by patient:

```
https://api.healthgorilla.com/fhir/Observation?  
category=https://www.healthgorilla.com/fhir/observation-category|func  
tional-status  
&patient=8063965523692105bfc35def
```

## 3. Create

HTTP Request	Method	Action
/Observation	POST	Creates Observation resource

The following attributes are required for Vital Sign:

Name	Details
subject	
related	Observation (code, value)

The following attributes are required for Social History:

Name	Details
subject	
code	
value	

The following attributes are required for Functional Status:

Name	Details
------	---------

subject	
code	
related	DeviceUseStatement

Example:

### 1. Add Vital Sign

```

<Observation xmlns="http://hl7.org/fhir">
  <meta>
    <profile
value="http://hl7.org/fhir/StructureDefinition/vitalspanel"/>
    </meta>
    <contained>
      <Observation xmlns="http://hl7.org/fhir">
        <id value="1"/>
        <status value="unknown"/>
        <category>
          <coding>
            <system
value="http://hl7.org/fhir/observation-category"/>
            <code value="vital-signs"/>
            <display value="Vital Signs"/>
          </coding>
        </category>
        <code>
          <coding>
            <system value="http://loinc.org"/>
            <code value="8310-5"/>
            <display value="Body temperature"/>
          </coding>
          <coding>
            <system value="http://snomed.info/sct"/>
            <code value="386725007"/>
            <display value="Body temperature"/>
          </coding>
        </code>
        <subject>
          <reference value="Patient/03db4352ee0190ad5e9e6748"/>
          <display value="Homer Simpson"/>
        </subject>
        <valueQuantity>
          <value value="36.5"/>
          <unit value="C"/>
          <system value="http://unitsofmeasure.org"/>

```

```

        <code value="Cel"/>
    </valueQuantity>
</Observation>
</contained>
<status value="unknown"/>
<category>
    <coding>
        <system
value="http://hl7.org/fhir/observation-category"/>
        <code value="vital-signs"/>
        <display value="Vital Signs"/>
    </coding>
</category>
<code>
    <coding>
        <system value="http://loinc.org"/>
        <code value="85353-1"/>
        <display value="Vital Signs Panel"/>
    </coding>
</code>
<subject>
    <reference value="Patient/03db4352ee0190ad5e9e6748"/>
    <display value="Homer Simpson"/>
</subject>
<effectiveDateTime value="2018-06-01"/>
<performer>
    <reference value="Practitioner/ceda4352e8f6405dd2e42f09"/>
    <display value="John Doe"/>
</performer>
<related>
    <type value="has-member"/>
    <target>
        <reference value="#1"/>
    </target>
</related>
</Observation>

```

## 2. Add Social History

```

<Observation xmlns="http://hl7.org/fhir">
    <status value="unknown"/>
    <category>
        <coding>
            <system
value="http://hl7.org/fhir/observation-category"/>
            <code value="social-history"/>

```

```

        <display value="Social History"/>
    </coding>
</category>
<code>
    <coding>
        <system value="http://snomed.info/sct"/>
        <code value="160573003"/>
        <display value="Alcohol Use"/>
    </coding>
    <text value="port wine, red"/>
</code>
<subject>
    <reference value="Patient/03db4352ee0190ad5e9e6748"/>
    <display value="Homer Simpson"/>
</subject>
<effectivePeriod>
    <start value="2017-12-06"/>
    <end value="2018-09-20T01:30:00+03:00"/>
</effectivePeriod>
<performer>
    <reference value="Practitioner/ceda4352e8f6405dd2e42f09"/>
    <display value="John Doe"/>
</performer>
<valueCodeableConcept>
    <coding>
        <display value="Only 10 old port wine"/>
    </coding>
</valueCodeableConcept>
    <comment value="This began just after the trip to Lisboa."/>
</Observation>

```

#### 4. Update

HTTP Request	Method	Action
/Observation/{ID}	PUT	Updates Observation resource

The following attributes are immutable for Vital Sign:

Name	Details
patient	
performer	

The following attributes are immutable for Social History:

Name	Details
patient	
performer	
code	

The following attributes are immutable for Functional Status:

Name	Details
patient	
performer	

## 5. Delete

HTTP Request	Method	Action
/Observation/{ID}	DELETE	Deletes Observation resource

## Procedure

The Procedure resource describes an action that is or was performed on a patient. This can be a physical intervention like an operation, or less invasive like counseling or hypnotherapy.

FHIR Resource: Procedure <https://www.hl7.org/fhir/procedure.html>

The following attributes are supported:

Name	Comment
<a href="#">status</a>	
<a href="#">category</a>	
<a href="#">code</a>	
<a href="#">subject</a>	Patient
<a href="#">performed</a>	
<a href="#">performer</a>	
<a href="#">location</a>	
<a href="#">usedReference</a>	Medication, Device
<a href="#">bodySite</a>	
<a href="#">reasonCode</a>	
<a href="#">reasonReference</a>	
<a href="#">outcome</a>	
<a href="#">followUp</a>	
<a href="#">partOf</a>	ProcedureRequest
<a href="#">note</a>	
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder">https://www.healthgorilla.com/fhir/StructureDefinition/record-recorder</a>	Recorder (Participant)
<a href="http://hl7.org/fhir/StructureDefinition/procedure-method">http://hl7.org/fhir/StructureDefinition/procedure-method</a>	CodeableConcept

## Endpoint

<https://api.healthgorilla.com/fhir/Procedure>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/Procedure/{ID}</code>	GET	Retrieve Procedure resource by ID

### 2. Search

HTTP Request	Method	Action
<code>/Procedure?[parameter={value}]</code>	GET	Retrieve Procedure resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier

Examples:

#### 1. Retrieve procedures by patient:

```
https://api.healthgorilla.com/fhir/Procedure?patient=8063965523692105bfc35def
```

### 3. Create

HTTP Request	Method	Action
<code>/Procedure</code>	POST	Creates Procedure resource

The following attributes are required:

Name	Details
<code>subject</code>	
<code>code</code>	

### 4. Update

HTTP Request	Method	Action
<code>/Procedure/{ID}</code>	PUT	Updates Procedure resource

The following attributes are immutable:

Name	Details
------	---------

patient	
recorder	
code	

## 5. Delete

HTTP Request	Method	Action
<code>/Procedure/{ID}</code>	DELETE	Deletes Procedure resource



## ProcedureRequest

The ProcedureRequest resource contains record of a request for diagnostic investigations, treatments, or operations to be performed.

FHIR Resource: ProcedureRequest <https://www.hl7.org/fhir/procedurerequest.html>

The following attributes are supported:

Name	Comment
<a href="#">identifier</a>	ACSN (Lab Reference ID)
<a href="#">requisition</a>	
<a href="#">status</a>	
<a href="#">intent</a>	
<a href="#">subject</a>	Patient
<a href="#">priority</a>	
<a href="#">category</a>	
<a href="#">performer</a>	Laboratory (Organization)
<a href="#">code</a>	
<a href="#">authoredOn</a>	
<a href="#">requester</a>	
<a href="#">reasonCode</a>	Diagnosis
<a href="#">supportingInfo</a>	AOE
<a href="#">specimen</a>	
<a href="#">note</a>	
<a href="http://hl7.org/fhir/StructureDefinition/procedurerequest-authorizedBy">http://hl7.org/fhir/StructureDefinition/procedurerequest-authorizedBy</a>	Author (Practitioner)

## Endpoint

<https://api.healthgorilla.com/fhir/ProcedureRequest>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/ProcedureRequest/{ID}</code>	GET	Retrieve ProcedureRequest by ID

### 2. Search

HTTP Request	Method	Action
<code>/Procedure?[parameter={value}]</code>	GET	Retrieve ProcedureRequest resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier
requisition	Token	Requisition Identifier
_lastUpdated	Date	

Examples:

#### 1. Retrieve ProcedureRequest resources by patient and date:

```
https://api.healthgorilla.com/fhir/Procedure?
patient=8063965523692105bfc35def
&_lastUpdated=2018-07
```

#### 2. Retrieve ProcedureRequest resources by order ID:

```
https://api.healthgorilla.com/fhir/Procedure?requisition=a5606c5b8998
949c08f49526
```

## RequestGroup

The RequestGroup resource represents a group of related requests that can be used to capture intended activities that have inter-dependencies. In HealthGorilla this resource used to describe a laboratory/radiology order.

FHIR Resource: RequestGroup <https://www.hl7.org/fhir/requestgroup.html>

The following profile is required:

```
https://healthgorilla.com/fhir/StructureDefinition/hg-order
```

The following attributes are supported:

Name	Comment
identifier	ACSN (Lab Reference ID), HG ID
requisition	
status	
intent	
subject	Patient
priority	
author	Practitioner
authoredOn	
action	ProcedureRequest
note	
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-authorizedBy">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-authorizedBy</a>	Author(Organization)
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-account">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-account</a>	Account
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-performer">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-performer</a>	Performer(Organization)
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-requester">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-requester</a>	Requester
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-location">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-location</a>	Room/Bed
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-reasoncode">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-reasoncode</a>	Diagnoses
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-specimen">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-specimen</a>	Specimen
<a href="https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-deliveryOptions">https://www.healthgorilla.com/fhir/StructureDefinition/requestgroup-deliveryOptions</a>	Delivery Options (location, method, electronic)
<a href="https://www.healthgorilla.com/fhir/StructureDe">https://www.healthgorilla.com/fhir/StructureDe</a>	Performer Location (Location)

inition/performer-location	
----------------------------	--

## Endpoint

<https://api.healthgorilla.com/fhir/RequestGroup>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/RequestGroup/{ID}</code>	GET	Retrieve RequestGroup by ID

### 2. Search

HTTP Request	Method	Action
<code>/RequestGroup?[parameter={value}]</code>	GET	Retrieve RequestGroup resources by the specified search criteria

Parameters:

Key	Type	Description
patient	Reference	Patient Identifier
identifier	Token	External Identifier
_lastUpdated	Date	

Examples:

#### 1. Retrieve Order by ID

<https://api.healthgorilla.com/fhir/RequestGroup/a5606c5b8998949c08f49526>

#### 2. Retrieve Order by requisition ID:

<https://api.healthgorilla.com/fhir/RequestGroup?identifier=https://www.healthgorilla.com/a5606c5b8998949c08f49526>

#### 3. Retrieve Orders by patient and date:

[https://api.healthgorilla.com/fhir/RequestGroup?patient=8063965523692105bfc35def&\\_lastUpdated=2018-07](https://api.healthgorilla.com/fhir/RequestGroup?patient=8063965523692105bfc35def&_lastUpdated=2018-07)

### 3. Create

Clients can submit new orders. Look for a separate documentation for that feature.

## DirectMessage

A DirectMessage resource is used to describe a message between persons and organizations in HealthGorilla Application. A message is a composite object that contains sender and recipients, subject, message body and attachments. Users may send messages to (and receive from) physicians, practice staff, patients, healthcare organizations.

This is a custom HealthGorilla resource, that is not defined by FHIR specification.

You can find detailed resource description in a separate specification.

The following attributes are supported:

Name	Comment
identifier	
sentOn	
sender	
sender.person	
sender.organization	
sender.telecom	
sender.display	
recipient	
recipient.person	
recipient.organization	
recipient.telecom	
recipient.display	
title	
content	
content.text	
content.html	
subject	
highPriority	

active	
attachment	
attachment.file	
attachment.reference	
owner	Practitioner
direction	IN / OUT

## Endpoint

<https://api.healthgorilla.com/fhir/DirectMessage>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
<code>/DirectMessage/{ID}</code>	GET	Retrieve DirectMessage by ID

### 2. Search

HTTP Request	Method	Action
<code>/DirectMessage?[parameter={value}]</code>	GET	Retrieve DirectMessage resources by the specified search criteria

Parameters:

Key	Type	Description
owner	Reference	Practitioner Identifier
direction	Token	
_lastUpdated	Date	

Examples:

1. Retrieve outgoing messages were sent by the given user since Sep'5:

```
https://api.healthgorilla.com/fhir/RequestGroup?
owner=8063965523692105bfc35def
&direction=OUT
&_lastUpdated=gt2018-09-05
```

## Subscriptions

The Subscription resource allows you to be notified about modifications of Health Gorilla resources in a most convenient way. You can use Subscription endpoint to create and manage subscriptions.

FHIR Resource: Subscription <https://www.hl7.org/fhir/subscription.html>

## Endpoint

<https://api.healthgorilla.com/fhir/Subscription>

## FHIR Operations

The following operations are currently supported:

### 1. Read

HTTP Request	Method	Action
/Subscription/{ID}	GET	Get Subscription by ID

### 2. Search

HTTP Request	Method	Action
/Subscription?[parameter={value}]	GET	

Parameters:

Key	Type	Description
status	Token	The current state of the subscription
type	Token	The type of channel for the sent notifications

### 3. Delete

HTTP Request	Method	Action
/Subscription/{ID}	DELETE	Deletes subscription

### 4. Create

HTTP Request	Method	Action
/Subscription	POST	Creates new subscription

Limited number of resources is supported, while some resources support filters that you can use in criteria:

Resource	Additional attributes you can use in criteria
DiagnosticReport	- patient - category

DocumentReference	- patient - class - type
RequestGroup	
ProcedureRequest	
DirectMessage	- owner - direction

## Examples

1. Retrieve Subscriptions by type:

<https://api.healthgorilla.com/fhir/Subscription?type=http://hl7.org/fhir/subscription-channel-type|rest-hook>

2. Create new subscription:

<https://api.healthgorilla.com/fhir/Subscription>

Headers:

Content-Type      application/json

Body:

```
{
  "resourceType": "Subscription",
  "status": "requested",
  "end": "2021-01-01T00:00:00Z",
  "reason": "Monitor new neonatal function",
  "criteria": "DiagnosticReport",
  "channel": {
    "type": "rest-hook",
    "endpoint": "http://localhost:8080/extsub/webhook",
    "payload": "application/fhir+json",
    "header": [
      "Authorization: Bearer secret-token-abc-123"
    ]
  }
}
```

## Subscription channel types

Health Gorilla supports rest-hook and web socket channels from the list

<https://www.hl7.org/fhir/valueset-subscription-channel-type.html>

## Rest hook

You can receive callback to the specified endpoint. And this is a highly recommended to add HTTP headers to the channel, that you will use to authorize the callback.

Only HTTPS endpoints can be used.



Supported payloads:

Payload	Description
No payload	See <a href="https://www.hl7.org/fhir/subscription.html#2.46.6.1">https://www.hl7.org/fhir/subscription.html#2.46.6.1</a> for details
application/fhir+xml	
application/fhir+json	
application/hg-event+json	Health Gorilla sends simple JSON event to client's endpoint with the resource type and resource's identifier that triggered the event. For example: <pre>{   "resource" : "RequestGroup",   "id" : "RequestGroup/d6fjlvj9gjb9djlaww8q" }</pre>

If you use no-payload mode, then you can query Health Gorilla FHIR server for the list of recently modified resources by appending “&\_lastUpdated=<date>” parameter to the search criteria.

## Web Socket

You can connect to secure Web Socket and receive runtime notifications.

To do so you need to use OAuth access token to authorize connection.

Web Socket endpoint is defined in the server's CapabilityStatement.

See <https://www.hl7.org/fhir/extension-capabilitystatement-websocket.html> for more information

# Error Handling

The API will always return sensible HTTP status code.

API errors are of 2 types:

- 400 series status codes for client issues
- 500 series status codes for server issues.

Common HTTP Status codes returned on errors:

- **400 Bad Request** - search could not be processed or failed basic FHIR validation rules
- **401 Not Authorized** - authorization is required for the interaction that was attempted
- **404 Not Found** - resource type not supported, or not a FHIR end-point
- **422 Unprocessable Entity** - resource was rejected by the server because it "violated applicable FHIR profiles or server business rules"
- **500 - Internal Server Error**

Along with the Status Code the API will return an OperationOutcome resource as response containing details about the error. Here is an example of OperationOutcome object:

```
{
  "resourceType": "OperationOutcome",
  "issue": [
    {
      "severity": "error",
      "code": "processing",
      "diagnostics": "Invalid request: The FHIR endpoint on this server does not know how to handle GET operation[Patient] with parameters [[wrong_parameter]]"
    }
  ]
}
```

## Appendix A

Resource samples

### DiagnosticReport

```
{
  "resourceType": "DiagnosticReport",
  "id": "f816835975440b570eed9683",
  "meta": {
    "lastUpdated": "2017-08-03T12:30:34.527+00:00"
  },
  "text": {
    "status": "generated",
```

```

    "div": "<div xmlns=\"http://www.w3.org/1999/xhtml\">Emr Test
Compass Lab Services Result 08/03/2017 12:30 PM</div>"
  },
  "contained": [
    {
      "resourceType": "ProcedureRequest",
      "id": "1",
      "status": "unknown",
      "intent": "order",
      "category": [
        {
          "coding": [
            {
              "system": "http://snomed.info/sct",
              "code": "103693007",
              "display": "Diagnostic procedure"
            }
          ]
        }
      ],
      "subject": {
        "reference": "Patient/6a1783599edf3ae52e51947a",
        "display": "Emr Test"
      },
      "requester": {
        "agent": {
          "reference":
"Practitioner/c50f6150e4b8c94c4ba3b23d",
          "display": "John Gorilla Doe, Ph. D. Psychology"
        }
      }
    }
  ],
  "identifier": [
    {
      "system": "https://www.healthgorilla.com",
      "value": "f816835975440b570eed9683"
    },
    {
      "type": {
        "coding": [
          {
            "system": "http://hl7.org/fhir/v2/0203",
            "code": "ACSN",
            "display": "Accession ID"
          }
        ]
      }
    }
  ]
}

```

```

        ],
        "text": "Accession Number"
    },
    "value": "770165"
}
],
"basedOn": [
    {
        "reference": "#1"
    }
],
"status": "final",
"category": {
    "coding": [
        {
            "system": "http://hl7.org/fhir/v2/0074",
            "code": "LAB",
            "display": "Laboratory"
        }
    ],
    "text": "Laboratory"
},
"code": {
    "text": "CMOR, PDF Report"
},
"subject": {
    "reference": "Patient/6a1783599edf3ae52e51947a",
    "display": "Emr Test"
},
"issued": "2017-08-02T17:49:10.000+00:00",
"performer": [
    {
        "actor": {
            "reference":
"Organization/f-f6c78159fc6d1019e302215f",
            "display": "Compass Lab Services"
        }
    }
],
"result": [
    {
        "reference": "Observation/f816835975440b570eed9683-2000",
        "display": "PDF Report"
    },
    {

```

```

        "reference":
"Observation/f816835975440b570eed9683-1304040",
        "display": "CMOR"
    }
]
}

```

## Observation (Observation Group)

```

{
  "resourceType": "Observation",
  "id": "f816835975440b570eed9683-1304040",
  "meta": {
    "lastUpdated": "2017-08-03T12:30:34.527+00:00"
  },
  "extension": [
    {
      "url":
"http://hl7.org/fhir/StructureDefinition/event-partOf",
      "valueReference": {
        "reference":
"DiagnosticReport/f816835975440b570eed9683",
        "display": "CMOR, PDF Report"
      }
    }
  ],
  "status": "unknown",
  "category": [
    {
      "coding": [
        {
          "system":
"http://hl7.org/fhir/observation-category",
          "code": "laboratory",
          "display": "Laboratory"
        }
      ]
    }
  ],
  "code": {
    "text": "CMOR"
  },
  "subject": {
    "display": "Emr Test"
  },
}

```

```

"effectiveDateTime": "2017-08-01T19:00:00+00:00",
"issued": "2017-08-02T17:49:10.000+00:00",
"performer": [
  {
    "reference": "Organization/f-f6c78159fc6d1019e302215f",
    "display": "Compass Lab Services"
  }
],
"related": [
  {
    "type": "has-member",
    "target": {
      "reference":
"Observation/f816835975440b570eed9683-1304040-1",
      "display": "Hydromorphone"
    }
  },
  {
    "type": "has-member",
    "target": {
      "reference":
"Observation/f816835975440b570eed9683-1304040-2",
      "display": "Hydromorphone, Quan"
    }
  },
  {
    "type": "has-member",
    "target": {
      "reference":
"Observation/f816835975440b570eed9683-1304040-3",
      "display": "Morphine"
    }
  },
  {
    "type": "has-member",
    "target": {
      "reference":
"Observation/f816835975440b570eed9683-1304040-4",
      "display": "Morphine, Quan"
    }
  }
]
}

```

## Observation (Observation Result)

```
{
  "resourceType": "Observation",
  "id": "f816835975440b570eed9683-1304040-1",
  "meta": {
    "lastUpdated": "2017-08-03T12:30:34.527+00:00"
  },
  "extension": [
    {
      "url":
"http://hl7.org/fhir/StructureDefinition/event-partOf",
      "valueReference": {
        "reference":
"Observation/f816835975440b570eed9683-1304040"
      }
    }
  ],
  "status": "final",
  "category": [
    {
      "coding": [
        {
          "system":
"http://hl7.org/fhir/observation-category",
          "code": "laboratory",
          "display": "Laboratory"
        }
      ]
    }
  ],
  "code": {
    "text": "Hydromorphone"
  },
  "subject": {
    "display": "Emr Test"
  },
  "effectiveDateTime": "2017-08-02T16:21:35+00:00",
  "issued": "2017-08-02T17:49:10.000+00:00",
  "valueString": "Negative",
  "interpretation": {
    "coding": [
      {
        "system": "http://hl7.org/fhir/v2/0078",
        "code": "N",
        "display": "Normal"
      }
    ]
  }
}
```

```
    }
  ]
},
"referenceRange": [
  {
    "text": "Cutoff level 50"
  }
]
}
```